

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

NNMT Antibody (Clone AT11G11)

单克隆抗体

产品货号: MA01576

第三版

免疫原: Recombinant human NNMT (1-264aa) purified from E. coli

亚型: Mouse IgG2a, k

克隆: Anti-human NNMT mAb, clone AT11G11, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human NNMT protein.

描述: N-methylation is one method by which drug and other xenobiotic compounds are metabolized by the liver. This gene encodes the protein responsible for this enzymatic activity which uses S-adenosyl methionine as the methyl donor.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF