



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

SET 7/9 Anti-Human Antibody (Clone s4E5)

单克隆抗体

产品货号: MA01577

第三版

免疫原: Recombinant human SET7/9 (1-366aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: s4E5

描述: SET7/9 is a histone methyltransferase(HMTase) that transfers methyl group to Lys of histone H3, in complex with S-adenosyl-L-methionine(AdoMet). The methylation of lysine residues of histones plays a critical role in the regulation of chromatin structure and gene expression.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF