

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

NPM Anti-Human Antibody (Clone 5E3)

单克隆抗体

产品货号: MA01583

第三版

免疫原: Recombinant human NPM (81-294aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: 5E3.

描述: Nucleophosmin (NPM; also known as B23, NO38, or numatrin) is a nucleolar phosphoprotein that binds the tumor suppressors p53 and p19Arf and is thought to be indispensable for ribogenesis, cell proliferation, and survival after DNA damage. Nucleophosmin is upregulated in many cancer cells.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF