

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

NQO2 Anti-Human Antibody (Clone AT7G7)

单克隆抗体

产品货号: MA01585

第三版

**免疫原:**Recombinant human NQO2 (1-231aa) purified from E. coli

**亚型:**Mouse IgG1 heavy chain and kappa light chain

**克隆:**AT7G7.

**描述:**Quinone oxidoreductase (NQO1 and NQO2) are cytosolic proteins that catalyze metabolic redaction of quinines and derivates.NQO2 is inhibited by flavones such as quercetin. Also benzo(a)pyrene is another known inhibitor of NQO2.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF