

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

NRG1(HRG- 3) Anti-Human Antibody (Clone k1G13)

单克隆抗体

产品货号: MA01586

第三版

**免疫原:**Recombinant human NRG1 isoform HRG-beta3 (20-241aa) purified from E. coli

**亚型:**Mouse IgG1 heavy chain and kappa light chain

**克隆:**k1G13

**描述:**Neuregulin1 (NRG1) is one of four proteins in the neuregulin family that act on the Epidermal growth factor receptor (EGFR) family of receptors. NRG1 was originally identified as a 44 kDa glycoprotein that interacts with the ERBB2 receptor tyrosine kinase to increase its phosphorylation on tyrosine residues.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF