

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Nsd3 Anti-Human Antibody (Clone 2E9)

单克隆抗体

产品货号: MA01588

第三版

免疫原: Recombinant human NSD3 (383-660aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: 2E9.

描述: This gene encoding NSD3 is related to the Wolf-Hirschhorn syndrome candidate - 1 gene and encodes a protein with PWWP (proline-tryptophan-tryptophan-proline) domains. NSD3 is Histone methyltransferase which preferentially methylates 'Lys-4' and 'Lys-27' of histone H3 and is amplified in several tumor cell lines and primary breast carcinomas.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF