

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-NSE (3-3-C) Mouse IgG MoAb

单克隆抗体

产品货号: MA01590

第三版

免疫原:NSE

亚型:IgG1

克隆:3-3-C

描述:Enolase is an enzyme in glycolysis system and the protein consists of dimer structure selected from three subunits  $\alpha$ ,  $\beta$  and  $\gamma$ . As  $\alpha$  and  $\beta$  Enolase exists such as in neuron cell and axis cylinder process, they are called “ Neuron specific Enolase (NSE) ”.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF