

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

OAT Anti-Human Antibody (AT23A2)

单克隆抗体

产品货号: MA01596

第三版

免疫原: Recombinant human OAT (33-439aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT23A2

描述: OAT is a 49-kDa nucleus-encoded protein imported into mitochondria to give the mature 48-kDa OAT polypeptide. It has been described in humans, animals, insects, plants and microorganisms. Especially OAT have sex-differential expression in the mouse kidney.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF