

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Occludin (1A8) Anti-Rat Mouse IgM MoAb

单克隆抗体

产品货号: MA01598

第三版

免疫原:Recombinant Rat Occludin

亚型:IgM

克隆:1A8

描述:Occludin is an approximately 60kDa-integral membrane protein localizing at Tight-Junction, and composed of short N-terminal cytosolic domain, four transmembrane domains and long C-terminal cytosolic domain.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF