

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

SNAP25 Anti-Human Antibody (Clone 4E11)

单克隆抗体

产品货号: MA01599

---

第三版

**免疫原:** Recombinant human SNAP25 (1-206aa) purified from E. coli

**亚型:** Mouse IgG1 heavy chain and kappa light chain

**克隆:** 4E+11

**描述:** The synaptosomal-associated protein (SNAP25) is an essential component of the core complex that mediates presynaptic vesicle trafficking. Thus, SNAP25 is directly involved in the release of neurotransmitters.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF