



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

OGG1 (7E2) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01601

第三版

免疫原:Recombinant Human OGG1

亚型:IgG1

克隆:700

描述:Oxidative damage to DNA has been proposed to have a role in cancer and ageing. Oxygen-free radicals formed during normal aerobic cellular metabolism attack bases in DNA, and 7, 8-dihydro-8-oxoguanine (8-oxoG) is one of the adducts formed.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF