

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

SNTA1 Antibody (Clone AT1E1)

单克隆抗体

产品货号: MA01602

第三版

免疫原: Recombinant human SNTA1 (1-505aa) purified from E. coli

亚型: IgG2b,k

克隆: Anti-human SNTA1 mAb, clone AT1E1, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human SNTA1 protein.

描述: SNTA1 (Syntrophin alpha 1) is 505 amino acid protein. SNTA1 is adapter protein that binds to and probably organizes the subcellular localization of a variety of membrane proteins. It may link various receptors to the actin cytoskeleton and the extracellular matrix via the dystrophin glycoprotein complex. SNTA1 plays an important role in synapse formation and in the organization of UTRN and acetylcholine receptors at the neuromuscular synapse. Defects of SNTA1 are the cause of long QT syndrome type 12 (LQT12).

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF