

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

SP-A (PE10) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01606

---

第三版

免疫原:Purified surfactant apoprotein A

亚型:IgG2b

克隆:PE10

描述:Lung surfactant is a lipoprotein complex which is synthesized and secreted into the alveoli of the lung by type II pneumocytes. Its major protein component, surfactant protein-A (SP-A), is a glycoprotein with a reduced denatured molecular mass of 35 kDa, and is a useful immunohistochemical marker for lung adenocarcinoma.

配方:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF