

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ORP150 (Oxygen Regulated Protein, p150) (2F07) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01613

第三版

免疫原:Recombinant Protein (GST-ORP150 (Leu 508 - Leu 999))

亚型:IgG2a

克隆:2F07

描述:Heat shock proteins (HSPs)/stress proteins are molecular chaperones that are induced by various environmental and physiological stimuli. The 150-kDa oxygen-regulated protein (ORP150), a new member of HSP family, functions as a molecular chaperone in the endoplasmic reticulum (1).

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF