



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Osteopontin (10A16) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01615

第三版

**免疫原:** Synthetic peptide for a part of Human Osteopontin

**亚型:** IgG1

**克隆:** 10A16

**描述:** Osteopontin (OPN) is a secreted glycoprotein that was originally isolated from bone. Its molecular weights have been reported in the range of 66 kDa to 44 kDa depending on glycosylation and phosphorylation.

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

**浓度:** Please see datasheet.

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF