
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Osteopontin (1B20) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01617

第三版

免疫原:Synthetic peptide in portion of C-terminus of Human Osteopontin

亚型:IgG1

克隆:1B20

描述:Osteopontin (OPN) is a secreted glycoprotein that was originally isolated from bone. Its molecular weights have been reported in the range of 66 kDa to 44 kDa depending on glycosylation and phosphorylation.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF