

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

SUMO2 Anti-Human Antibody (AT10F1)

单克隆抗体

产品货号: MA01620

第三版

**免疫原:**Recombinant human SUMO2 (1-93aa) purified from E. coli

**亚型:**Mouse IgG2b heavy chain and kappa light chain

**克隆:**AT10F1.

**描述:**SUMO (small ubiquitin-related modifier) is the best-characterized member of a growing family of ubiquitin-related proteins. 8 – 11 kDa that covalently modify various intracellular proteins. It resembles ubiquitin in its structure, its ability to be ligated to other proteins.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF