

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti- α -Synuclein (61-95 aa) Antibody (Clone 5C2)

单克隆抗体

产品货号: MA01622

第三版

免疫原: Recombinant human α -synuclein (61-95) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 5C2

描述: α -synuclein (amino acids 1-140), an acidic neuronal protein of 140 amino acids, is extremely heat-resistant and is natively unfolded with an extended structure primarily composed of random coils. α -synuclein has been suggested to be implicated in the pathogenesis of Parkinson's disease and related neurodegenerative disorders, and more recently, to be an important regulatory component of vesicular transport in neuronal cells.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF