



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Osteopontin N-Half (34E3) Mouse IgG MoAb

单克隆抗体

产品货号: MA01623

第三版

免疫原: Synthetic peptide for a part of Human Osteopontin (SVVYGLR)

亚型: IgG1

克隆: 34000

描述: Osteopontin (OPN) is a secreted glycoprotein that was originally isolated from bone. Its molecular weights have been reported in the range of 66 kDa to 44 kDa depending on glycosylation and phosphorylation.

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF