

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

p53 (Bp53-12) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01631

---

第三版

免疫原:Recombinant Human p53

亚型:IgG2a

克隆:Bp53-12

描述:Biochemically, p53 acts as a transcription factor that works in the form of tetramer (it is formed by the interaction of dimers that combined in C-end area), and binds to specific base sequence and then activate the transcription.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF