

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PAI1 Anti-Human Antibody (Clone AT1C5)

单克隆抗体

产品货号: MA01634

第三版

**免疫原:** Recombinant human PAI1 (24-402aa) purified from E. coli.

**亚型:** Mouse IgG1 heavy chain and kappa light chain

**克隆:** Anti-human PAI1 mAb, clone AT1C5, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PAI1 protein.

**描述:** Plasminogen activator inhibitor-1 (PAI-1), a member of the serine protease inhibitor (serpin) superfamily. PAI-1 is the principal inhibitor of tissue plasminogen activator and urokinase, the activators of plasminogen and hence fibrinolysis. PAI-1 is mainly produced by the endothelium but is also secreted by other tissue types, such as adipose tissue.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol

**浓度:** 1mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF