

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PAIP2 Anti-Human Antibody (Clone AT1C12)

单克隆抗体

产品货号: MA01636

第三版

免疫原:Recombinant human PAIP2 (1-127aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:AT1C12

描述:PAIP2 (Poly(A)-binding protein-interacting protein 2) is a 127 amino acid protein. PAIP binds to the polyadenylate-binding protein (PABP). There are two PAIP proteins called PAIP1 and PAIP2. PAIP2 acts as a repressor in the regulation of translation initiation of poly(A)-containing mRNAs. and Paip1 stimulates translation. PAIP2 contains two binding sites for PABP, one encompassing a 16-amino-acid stretch located in the C terminus and a second encompassing a larger central region.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF