

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Park7/DJ-1 Anti-Human Antibody (Clone 1B11)

单克隆抗体

产品货号: MA01640

第三版

**免疫原:**Recombinant human Park7/DJ-1 (1-189 aa) purified from E. coli

**亚型:**Mouse IgG1 heavy chain and kappa light chain

**克隆:**1B11

**描述:**Park7/DJ-1, which is a small 189 amino acid protein, is a ubiquitously expressed protein involved in various cellular processes including spermatogenesis and fertilization, cancer, RNA-binding, androgen-receptor signaling and oxidative stress.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF