
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Parkin (5A1) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01647

第三版

免疫原: Synthetic peptide of part of C terminal of human parkin (AYRVDERRAAEQARWEAA)

亚型: IgG1

克隆: 5A1

描述: Parkinson's Disease (PD) is a relatively common neurodegenerative disorder, which is characterized by the loss of midbrain dopamine (DA) neurons and the presence of Lewy bodies, proteinaceous cytoplasmic inclusions that are abundantly enriched in ubiquitin.

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF