

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PCBP1 Antibody (Clone AT2A10)

单克隆抗体

产品货号: MA01653

第三版

免疫原: Recombinant human PCBP1 (1-163aa) purified from E. coli

亚型: IgG2a ,k

克隆: Anti-human PCBP1 mAb, clone AT2A10, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PCBP1 protein.

描述: PCBP-1 and PCBP-2 have paralogues (PCBP3 and PCBP4) which are thought to have arisen as a result of duplication events of entire genes. PCBP1 constitutes a set of polypeptides that contribute to mRNA transcription, pre-mRNA processing as well as mature mRNA transport to the cytoplasm and translation. It has also been implicated in translational control of the 15-lipoxygenase mRNA, human Papillomavirus type 16 L2 mRNA, and hepatitis A virus RNA.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF