

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Pcna Anti-Human Antibody (Clone K1A3)

单克隆抗体

产品货号: MA01655

第三版

**免疫原:**Recombinant human PCNA (1-261aa) purified from E. coli

**亚型:**Mouse IgG3 heavy chain and kappa light chain

**克隆:**k1A3

**描述:**PCNA (Proliferating Cell Nuclear Antigen) is found in the nucleus and is a cofactor of DNA polymerase delta. This protein is associated with DNA synthesis and repair. The encoded protein acts as a homotrimer and helps increase the processivity of leading strand synthesis during DNA replication.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF