

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Pdcd1 Anti-Human Antibody (Clone 4F12)

单克隆抗体

产品货号: MA01656

第三版

免疫原:Recombinant human PDCCD1 (21-167aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:4F12

描述:PDCCD1 is an immunoreceptor that belongs to the immunoglobulin (Ig) superfamily, also known as Programmed Cell Death 1. It is expressed on activated T cell, B cells, myeloid cells and a subset of thymocytes.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF