
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PDCD12 Anti-Human Antibody (Clone 3G4)

单克隆抗体

产品货号: MA01657

第三版

免疫原: Recombinant human PDCD12 (254-362aa) purified from E. coli

亚型: Mouse IgG2a heavy chain and kappa light chain

克隆: 3G4

描述: PDCD12 (Programmed cell death 12), also referred to as AVEN, represents a new class of cell death regulator. PDCD12 is a conserved intracellular membrane protein that blocks apoptosis impeded by Apaf-1 mediated caspase activation.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF