

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Pdcd4 Anti-Human Antibody (Clone K4C1)

单克隆抗体

产品货号: MA01659

第三版

免疫原:Recombinant human PDCCD4 (1-469aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:K4C1

描述:Programmed cell death 4 (PDCCD4) encodes a tumor suppressor protein whose expression is lost in progressed carcinomas of lung, breast, colon, and prostate. The expression of PDCCD4 gene is strongly induced during apoptosis in a number of cell types. It is modulated by cytokines in natural killer and T cells.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF