

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PDK1 Anti-Human Antibody (Clone AT2D6)

单克隆抗体

产品货号: MA01667

第三版

免疫原: Recombinant human PDK1 (29-436aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT2D6

描述: PDK1 (pyruvate dehydrogenase kinase isoform 1) is involved in the regulation of enzymatic activity of mammalian pyruvate dehydrogenase (PDH) that is a part of a mitochondrial multienzyme complex to catalyze the oxidative decarboxylation of pyruvate and is one of the major enzymes responsible for the regulation of homeostasis of carbohydrate fuels in mammal

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF