



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Pdpn Anti-Human Antibody(Clone 5E2)

单克隆抗体

产品货号: MA01668

第三版

免疫原:Recombinant human PDPN (1-206aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:500

描述:Podoplanin(PDPN) is a small mucin-like transmembrane protein, widely expressed in various specialized cell types throughout the body. This is a type-I integral membrane glycoprotein with diverse distribution in human tissues. The physiological function of this protein may be related to its mucin-type character.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF