



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PDPN Antibody (Clone AT56F7)

单克隆抗体

产品货号: MA01670

第三版

免疫原: Recombinant human PDPN (99-207aa) purified from E. coli

亚型: Mouse IgG1, k

克隆: Anti-human PDPN mAb, clone AT56F7, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PDPN protein.

描述: Podoplanin (PDPN) is a small mucin-like transmembrane protein, widely expressed in various specialized cell types throughout the body. This is a type-I integral membrane glycoprotein with diverse distribution in human tissues. The physiological function of this protein may be related to its mucin-type character. The homologous protein in other species has been described as a differentiation antigen and influenza-virus receptor.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF