



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PEDF Anti-Human Antibody (Clone 1C4)

单克隆抗体

产品货号: MA01674

第三版

免疫原: Recombinant human PEDF (19-418 aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 1C4

描述: Pigment epithelium-derived factor (PEDF) is a 50kDa glycoprotein and member of the serpin superfamily without any activity of serine protease. PEDF is an angiogenesis inhibitor with neurotrophic properties, balances angiogenesis in the eye and blocks tumor progression.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF