



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PGP9.5 Anti-Human Antibody (Clone 3D9)

单克隆抗体

产品货号: MA01691

第三版

**免疫原:** Recombinant human PGP9.5 (1-223aa) purified from E. coli

**亚型:** Mouse IgG2a heavy chain and kappa light chain

**克隆:** 3D9

**描述:** PGP9.5/UCH-L1 is a member of a gene family whose products hydrolyze small C-terminal adducts of ubiquitin to generate the ubiquitin monomer. PGP9.5 is a component of the ubiquitin system, which has a fundamental role in regulating various biological activities. PGP9.5 gene encodes two opposing enzymatic activities that affect alpha-synuclein degradation and Parkinson's disease susceptibility.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF