

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Pin1 Anti-Human Antibody (Clone 3G8)

单克隆抗体

产品货号: MA01693

第三版

**免疫原:**Recombinant human Pin1 (1-163aa) purified from E. coli

**亚型:**Mouse IgG1 heavy chain and kappa light chain

**克隆:**3G8

**描述:**Human Pin1 is a peptidyl-prolyl cis/trans isomerase (PPIase) that interacts with NIMA and essential for cell cycle regulation. Pin1 is nuclear PPIase containing a WW protein interaction domain, and is structurally and functionally related to Ess1/Ptf1, an essential protein in budding yeast.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF