

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PKLR Anti-Human Antibody (Clone AT1E3)

单克隆抗体

产品货号: MA01695

第三版

免疫原: Recombinant human PKLR (47-574aa) purified from E. coli

亚型: IgG1 heavy chain and kappa light chain.

克隆: AT1E3

描述: PKLR is a pyruvate kinase that catalyzes the transphosphorylation of phosphoenolpyruvate into pyruvate and ATP, which is the rate-limiting step of glycolysis. The PKLR gene encodes the L- and R-type isoenzymes through alternative splicing events under the control of different promoters.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF