



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PNMT Anti-Human Antibody (Clone AT1C11)

单克隆抗体

产品货号: MA01705

第三版

免疫原:Recombinant human PNMT (1-282aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and light chain

克隆:AT1C11

描述:Phenylethanolamine N-methyltransferase (PNMT) catalyzes the synthesis of epinephrine and plays an important role in regulating cardiovascular function. Expression of PNMT is regulated by hormonal and neural stimuli because its promoter contains sequences responding to cholinergic and depolarization stimuli.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在-20 °C 或者-70 °C. 避免反复冻融.

应用:ELISA, WB, ICC/IF