

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PNPO Anti-Human Antibody (Clone AT2C7)

单克隆抗体

产品货号: MA01706

第三版

免疫原: Recombinant human PNPO (57-261aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT2C7

描述: PNPO (pyridoxamine 5'-phosphate oxidase) is a 261 amino acid protein belonging to the pyridoxamine 5'-phosphate oxidase family. It is the rate-limiting enzyme in vitamin B6 synthesis. Vitamin B6 (Pyridoxal 5-prime-phosphate or PLP) is vital for normal cellular function, and some cancer cells have notable differences in vitamin B6 metabolism compared to their normal counterparts.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF