

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PP2C Anti-Human Antibody (Clone p6c7)

单克隆抗体

产品货号: MA01712

第三版

免疫原: Recombinant human PP2C (1-382aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: p6c7

描述: Protein phosphatase 2C (PP2C) is a Mn^{2+} - or Mg^{2+} -dependent protein serine/threonine phosphatase that inhibits the human stress-responsive p38 and JNK MAPK pathways and regulates cellular stress response in eukaryotes.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF