



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PPIC Anti-Human Antibody (Clone AT3C6)

单克隆抗体

产品货号: MA01716

第三版

**免疫原:** Recombinant human PPIC (39-197aa) protein purified from E. coli

**亚型:** IgG2b, k

**克隆:** AT3C6

**描述:** PPIC is a member of the peptidyl-prolyl cis-trans isomerase (PPIase) family. PPIases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and accelerate the folding of proteins. Similar to other PPIases, this protein can bind immunosuppressant cyclosporin A.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF