



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Ppm1B Anti-Human Antibody (Clone K1B1)

单克隆抗体

产品货号: MA01722

第三版

免疫原: Recombinant human PPM1B (1-192aa) purified from E. coli

亚型: Mouse IgG2a heavy chain and kappa light chain

克隆: K1B1

描述: Protein phosphatase 1B (PPM1B) is a member of the PP2C family of serine/threonine protein phosphatases. PP2C family members are known to be negative regulators of cell stress response pathways.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF