

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PPM1G Anti-Human Antibody (Clone k1G6)

单克隆抗体

产品货号: MA01725

第三版

免疫原: Recombinant human PPM1G (317-546aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: k1G6

描述: PPM1G (Protein phosphatase 1G) is a member of the PP2C family of Serine/threonine protein phosphatases. PP2C family members are known to be negative regulators of cell stress response pathways. The human Ser/Thr phosphatase PPM1G promotes spliceosome assembly.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF