

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Ppp1Ca Anti-Human Antibody (Clone 4G3)

单克隆抗体

产品货号: MA01727

第三版

**免疫原:**Recombinant human PPP1CA (30-299aa) purified from E. coli

**亚型:**Mouse IgG2a heavy chain and kappa light chain

**克隆:**4G3

**描述:**Protein phosphatase 1, catalytic subunit, alpha isoform, also known as PPP1CA, is one of the three catalytic subunits of protein phosphatase 1. The protein is highly conserved in eukaryotic cells and play a role in dephosphorylating substrates such as the postsynaptic density-associated Ca<sup>2+</sup>/calmodulin dependent protein kinase .

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF