



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Ppp1R14A Anti-Human Antibody (Clone 4H10)

单克隆抗体

产品货号: MA01729

第三版

免疫原: Recombinant human PPP1R14A (1-147aa) purified from E. coli

亚型: Mouse IgG2a heavy chain and kappa light chain

克隆: 4H10

描述: PPP1R14A (Protein phosphatase 1 regulatory subunit 14A) is a phosphorylation-dependent inhibitory protein for smooth muscle myosin phosphate. Myosin phosphatase can reverse MYL (myosin light chain) phosphorylation to induce a state of relaxation.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF