



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Presenilin-1 (17C2) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01738

第三版

免疫原:N-terminal fragment of recombinant Human Presenilin-1 (E. coli)

亚型:IgG1

克隆:17C2

描述:The presenilin 1 (PS1) gene was identified as the gene that harbors mutations that cause this type of FAD (termed AD3). PS1 is a novel 467-aa long integral membrane protein localized to intracellular membranous compartments, especially the endoplasmic reticulum, that spans the membrane 7-9 times.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF