

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Rab5A Anti-Human Antibody (Clone 3A4)

单克隆抗体

产品货号: MA01753

第三版

免疫原:Recombinant human Rab5a (1-215aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:3A4

描述:Rab5a is a member of the RAS oncogene family. Rab5a is a small GTPase localized to the plasma membrane and early endosomes. It is a crucial regulatory component of the endocytic pathway. Activation of Rab5a is mediated by GDP-GTP exchange factors (GEFs) that generate the Rab5-GTP complex.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF