



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Galectin 3 Monoclonal Antibody

单克隆抗体

产品货号: MA01758

第三版

免疫原:Synthetic Peptide

亚型:IgG

克隆:1A2

描述:Galactose-specific lectin which binds IgE. May mediate with the alpha-3, beta-1 integrin the stimulation by CSPG4 of endothelial cells migration. Together with DMBT1, required for terminal differentiation of columnar epithelial cells during early embryogenesis.

配方:

浓度:1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF