



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CDX2 Monoclonal Antibody

单克隆抗体

产品货号: MA01771

第三版

免疫原:Synthetic Peptide

亚型:IgG

克隆:1A5

描述:This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded protein is a major regulator of intestine-specific genes involved in cell growth and differentiation. This protein also plays a role in early embryonic development of the intestinal tract. Aberrant expression of this gene is associated with intestinal inflammation and tumorigenesis.

配方:

浓度:1mg/mL

保存: +4 °C 保存 (1-2 周). 长期保存在-20 °C 或者-70 °C. 避免反复冻融.

应用:ELISA, WB, ICC/IF