

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD66e Monoclonal Antibody

单克隆抗体

产品货号: MA01775

第三版

免疫原:Synthetic Peptide

亚型:IgG

克隆:1A2

描述:Carcinoembryonic antigen-related cell adhesion molecule 5 (CEACAM5) also known as CD66e (Cluster of Differentiation 66e), is a human gene. Members of the CEACAM subfamily, including CEACAM5, belong to the CEA gene family. CD66E is a well known tumor marker and a heavily glycosylated GPI-linked cell surface molecule.

配方:

浓度:1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF