



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CASP9 Monoclonal Antibody

单克隆抗体

产品货号: MA01786

第三版

免疫原:Synthetic Peptide

亚型:IgG

克隆:3I3

描述:Involved in the activation cascade of caspases responsible for apoptosis execution. Binding of caspase-9 to Apaf-1 leads to activation of the protease which then cleaves and activates caspase-3. Proteolytically cleaves poly(ADP-ribose) polymerase (PARP). Isoform 2 lacks activity is an dominant-negative inhibitor of caspase-9.

配方:

浓度:1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF